

Our strengths

Fast analysis: 1hr characterisation, sorting up to 100,000 cells per second

Accurate and reliable results: wide resolution range from 0.2 to 30 µm

Multiparametric sorting and analysis
Sorter solely dedicated to microorganisms

Development of custom methodologies and processes: reporter choice, sorting processes, ...

Dedicated specialised scientist

Strain engineering Analytical Biotransformation & culture Ethics

An integrated solution tailored to the needs, to accelerate the development of industrial biotechnology

Our service

- Individual quantitative and qualitative characterisation of microorganisms in suspension
- > Very high-speed cell sorting
- Prokaryotic or eukaryotic cells in suspension
- Direct analysis of the intrinsic optical properties and fluorescence induced
- Direct and standardised numeration with viability and vitality staining
- Selection of homogeneous populations

Deliverables

- > Enumeration
- Physiological data; intracellular pH, oxidative stress, membrane potential and membrane fluidity
- > Validation of a gene alteration
- > Isolated cellular populations in suitable culture tubes (up to 6)
- > Isolated clones in plates of 96, 384 or 1536 wells with one cell per well
- > Analysis reports (cytograms and statistics)


Equipment

- > Cell sorter MoFlo Astrios EQ, Beckman Coulter (5 lasers, jet-in-air, dual forward scatter technology, fully integrated in BSLII)
- > Cytometer analyser MACSQuant VYB, Miltenyi Biotec (3 lasers, absolute counts performance, plate autosampler)
- > FlowJo v10.1 (Tree Star, Inc): reference software for advanced flow cytometry data analysis



Some achievements

- Population assessment in microbial consortia
- **Detection and physical isolation** of viable but non-culturable populations
- > Detailed analysis of microbial catalysts for process improvement
- > **High throughput screening** of cellular banks and selection of mutants of interest



Additional services

- Transcriptomic and metabolomic dynamic studies
- Scale-up: microbial culture, enzymatic reactor

Contact

Manager: **Delphine Lestrade** E-mail: **delphine.lestrade@inra.fr**

. Tel.: **+33 (0)5 61 28 55 82**

Customised estimate on request

